













At a glance Disposable lenses


DISPOSABLES									
Minus BVP Range	PRODUCT NAME	BOZR mm	TD mm	H ₂ O %	CT mm @-3.00	Dk	CYL	AXES/ OTHER	Plus BVP Range


DAILIES										
-10	-6	DAILIES AquaComfort Plus® (non-ionic) 	8.70	14.00	69	0.10	26			+6
-10	-6	FOCUS® DAILIES® with AquaComfort® (non-ionic) 	8.60	13.80	69	0.10	26			+6
-10	-6	FOCUS® DAILIES® BASIC (non-ionic) 	8.60	13.80	69	0.10	26			+6
-8	-6	FOCUS® DAILIES® TORIC with AquaComfort® (non-ionic) 	8.60	14.20	69	0.10	26	-0.75 -1.50	90°&180°	+4
	-6	FOCUS® DAILIES® PROGRESSIVES with AquaComfort® (non-ionic) 	8.60	13.80	69	0.11	26		ADD UP TO +3.00	+5
	-6	FRESHLOOK® ONE-DAY (non-ionic) 	8.60	13.80	69	0.10	26			






SILICONE HYDROGEL										
-10	-8	AIR OPTIX® (non-ionic) 	8.60	14.20	33	0.08	110			+6
	-20	AIR OPTIX Individual™ (non-ionic) 	7.4, 7.7, 8.0, 8.3 7.8, 8.1, 8.4, 8.7, 9.0 8.0, 8.3, 8.6, 8.9, 9.2	13.2 14.0 14.8	32	0.07	82			+20
	-6	AIR OPTIX® for ASTIGMATISM (non-ionic) 	8.70	14.50	33	0.102	110	-0.75 -1.25	10° ATC	
	-6	Note: not all parameters are available at launch	8.70	14.50	33	0.102	110	-0.75 -1.25 -1.75 -2.25	10° ATC	
	-9	Parameters to follow	8.70	14.50	33	0.102	110	-1.25 -1.75	180°±20 and 90°±20 (10)	
	-9		8.70	14.50	33	0.102	110	-0.75 -1.25 -1.75		+6
-10	-8	AIR OPTIX® NIGHT&DAY® (non-ionic) 	8.40 8.60	13.80	24	0.08	140			+6

MONTHLY REPLACEMENT										
-15	-8	FOCUS® VISITINT® 	8.60 8.90	14.00	55	0.10	20			+6
-10	-6	PRECISION UV® (non-ionic) 	8.40	14.40	74	0.14	38			+5
-16	-6		8.70	14.40	74	0.14			+5	
	-6		9.10	14.40	74	0.14			+10	

DISPOSABLES									
Minus BVP Range	PRODUCT NAME	BOZR mm	TD mm	H ₂ O %	CT mm @-3.00	Dk	CYL	AXES/ OTHER	Plus BVP Range

TORIC									
-6	FOCUS® TORIC 	8.90 9.20	14.50	55	0.146	20	-1.00 -1.75	Minus BVP 10° ATC	
-6		8.90	14.50	55	0.146	20	-2.50	180°±20 and 90°±20 (10)	
		8.90	14.50	55	0.146	20	-1.00 -1.75	Plus BVP 180°±20 and 90°±20 (10)	+4

MULTIFOCAL									
-7	FOCUS PROGRESSIVES® 	8.60	14.00	55	0.10	16		ADD UP TO +3.00	+6

COLOURS*									
-8	FRESHLOOK® COLORS* 	8.60	14.50	55	0.08	16			+6
-8	FRESHLOOK COLORBLENDS®** 	8.60	14.50	55	0.08	16			+6
-8	FRESHLOOK DIMENSIONS®** 	8.60	14.50	55	0.08	16			+6
-8	FRESHLOOK RADIANCE®** 	8.60	14.50	55	0.08	16			+6
-8	FOCUS® SOFTCOLORS* 	8.60 8.90	14.00	55	0.10	20			+6

NEW
FreshLook® range available:
- 6.00D to +2.00D
in 0.25D steps
- 6.00D to -8.00D
in 0.50D steps
+ 2.00D to +6.00D
in 0.50D steps

NOVELTY*									
WILDEYES® 2 GO discontinued 		8.60	14.50	55	0.08	16		PLANO ONLY	

Coding key for powers ■ 0.25 steps ■ 0.50 steps

*For colour availability and illustrations please see pages 12 and 13.